

AED MAINTENANCE CHECKLIST

Date								
AED Unit Clean, no dirt or contamination; no damage present								
Supplies Available: a. 2 sets defibrillation pads, within expiration date, undamaged b. 1 set of pediatric defibrillation pads within expiration date, undamaged c. Ancillary supplies: <ul style="list-style-type: none"> • Towel • Scissors • Razor • CPR mask • Gloves d. Store unopened battery within install before date e. PC data card undamaged*								
Status Indicator:								
Inspected by:								
Remarks, Problems, Corrective Actions:								

Prepared by Emergency Skills, Inc.
 Daily, check AED status indicator.
 Monthly, perform thorough check, Date and initial upon completion.
 Please make copies of this form as needed. Store completed forms in a safe place.
 To order supplies, call the Dept of Education - Sherry Fennessy (212) 374-2320 or Donovan Lauther at (212) 374-6757.